## Calibrating the Main Temperature Control

We recommend that you calibrate your unit once it has been installed in its working environment and has been stable at the set point for 24 hours. To calibrate the unit, perform the following steps:

1. Place a certified reference thermometer in the chamber by placing it either directly inside or through the access tube at the top left of the unit. Ensure that the thermometer is not touching any shelving. If you place the thermometer directly inside the chamber, taping the thermometer to a Petri dish will raise it off the shelf and keep the scale in view.
2. Allow the temperature to stabilize again until the thermometer reads a constant value for one hour.
3. Compare the digital display with the reference thermometer.
4. If there is an unacceptable difference, put the display into calibration mode by pressing both the Up and Down arrow pads at the same time for approximately five (5) seconds until the display blinks off and on.
5. While blinking, the display can be calibrated by pressing the Up or Down arrow pads until the display reads the correct value.
6. Allow the incubator temperature to stabilize again, and recalibrate if necessary.

## Setting the Overtemperature Thermostat

To set the Overtemperature Thermostat, perform the following steps.

1. Once the incubator is stable at the desired set point, turn the Overtemperature Thermostat counterclockwise until the Overtemperature indicator light turns on.
2. Turn the Overtemperature Thermostat clockwise just the safety indicator light turns off.
3. Continue to turn the Overtemperature Thermostat clockwise an additional two small divisions on its scale past the point where the indicator light went out. This will set the Safety Thermostat at approximately 1°C above the Main Temperature set point.

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